




PHOENIX FLOW SYSTEMS

Flow Cytometer Sheath Fluid Concentrate



VS



2 LBS
fluid concentrate
of flow cytometer
sheath solution

**CALL
FOR A
FREE
SAMPLE**
1.800.886.FLOW

THE BENEFITS OF CHEAP SHEATH

- ⇒ **save on shipping** 2 lbs vs. 44 lbs
- ⇒ **save on storage space** ships in a 1 qt. bottle
- ⇒ **environmentally friendly** less land fill waste
- ⇒ **easy to formulate** just add water
- ⇒ **less expensive** from start to finish

6790 Top Gun St. #1 San Diego, CA 92121 call 1.800.886.FLOW
or 858.453.5095 fax 858.453.2117 email info@phoenixflow.com
or visit us online at phoenixflow.com

Cubitainer is a registered trademark of Hedwin Corp.

PHOENIX FLOW SYSTEMS

Cheap Sheath is a concentrated pH balanced formulation of Phosphate Buffered Saline which has been pre-filtered and is ready to use in any flow cytometer as your sheath solution upon dilution with distilled water. Each bottle of Cheap Sheath contains enough concentrate to prepare 20 liters of flow cytometer sheath solution.

Formulation is simple for flow cytometers that do not use 20 liter sheath containers.

For any flow cytometer with a 3.5 liter sheath tank, add 100 mls Cheap Sheath to 2.9 liters of distilled water in the sheath tank.